

Symptom-Rhythm Correlation Patterns in Patients Undergoing Ambulatory ECG Monitoring: Analysis of over 1 Million Patients

Anthony J. Battisti PhD¹, Relana Pinkerton PhD¹, Vladimir Fokin PhD¹, Brent Wright RN, DrPHc¹, Mintu P. Turakhia MD, MAS^{1,2}

¹iRhythm Technologies, San Francisco, CA; ²Stanford University, Stanford, CA

Background

- Symptoms are the most common indication for ambulatory cardiac monitoring.
- However, symptom-rhythm correlation (SRC) has not been well described, particularly across the spectrum of different rhythms.


Objective

- We assessed SRC in patients undergoing patch-based long-term continuous ambulatory ECG monitoring (LTCM) up to 14 days.

Methods

- Retrospective cohort study of Zio[®] monitor or XT LTCM (iRhythm Technologies, San Francisco, CA) devices prescribed in the US between June 2023 and July 2024 (**Figure 1**).
- Devices included a patient-activated button to document symptoms during the continuous recording.
- Inclusions: Patients ≥ 18 years; Devices worn for >7 to 14 days.
- Exclusions: 100% AF.
- Symptoms documented within a ± 45 second window of an arrhythmia were considered rhythm-correlated.
- For each rhythm type, a patient was classified as having SRC if there was ≥ 1 rhythm-correlated button press during the wear period, regardless of the number of button presses that were not temporally associated with an arrhythmia.
- Data were further stratified by observed symptom frequency, defined as the number of button presses/day:
 - Daily Symptoms (≥ 1 /day), or
 - Non-Daily Symptoms (< 1 /day)
- ECG data was analyzed via a deep-learned AI algorithm and confirmed by certified cardiographic technicians.

Results



Cohort Description

- Cohort size: 1,100,337; age 61.4 ± 17.6 years; 55.6% female.
- 20.8% (229,005) had Daily Symptoms.
- 79.2% (871,332) had Non-Daily Symptoms, including 284,784 patients with no documented symptoms.
- 78.5% had one or more arrhythmias (AF, SVT, VT, VF, Pause or 2nd Degree or Complete Heart Block):
 - AF (7.2%), SVT (71.3%), VT (25.5%), VF (0.01%), Pause (3.2%) and 2nd Degree or Complete Heart Block (1.5%).
- Most patients had atrial (99.6%) or ventricular (98.8%) ectopic beats.

SRC Findings

- There was wide variation in SRC by arrhythmia type and symptom frequency (**Figure 2**).
- SRC was highest for AF, and isolated ectopic beats (SVE and VE).

- Among patients with AF, SRC was 43.1%, including:
 - 83.8% SRC for patients with Daily Symptoms, and
 - 37.6% SRC for patients with Non-Daily Symptoms
- SRC was 34.9% for isolated SVE, including:
 - 67.2% SRC for patients with Daily Symptoms, and
 - 26.4% SRC for patients with Non-Daily Symptoms
- SRC was 32.5% for isolated VE, including:
 - 64.4% SRC for patients with Daily Symptoms, and
 - 24.2% SRC for patients with Non-Daily Symptoms
- SRC varied, but was low ($< 15\%$) for patients with other arrhythmias or PAC/PVC morphologies.
- SRC was higher for the Daily Symptom group vs. the Non-Daily Symptom group across all rhythms.
- Among patients with any arrhythmia detected, SRC was low (18.3%), including:
 - 43.1% SRC for patients with Daily Symptoms, and
 - 13.2% SRC for patients with Non-Daily Symptoms

Limitations

- SRC was defined as requiring only one temporally associated arrhythmia. This could have been coincidental in patients who had a high volume of non-correlated button presses.
- Symptoms may have been present but not reported by button press in some cases.


Conclusions

- This is the largest study characterizing symptom-rhythm correlation by rhythm type in ambulatory monitoring.
- The results demonstrate that **AF and ectopy** were the rhythms most correlated with patient symptoms.
- **Overall symptom-rhythm correlation was low**, but higher for patients with Daily Symptoms than Non-Daily Symptoms.

Disclosures

- AJ Battisti, R Pinkerton, V Fokin, and B Wright are employees of and have received equity from iRhythm Technologies, Inc.
- Dr. Turakhia has received equity from iRhythm, Connect America, Evidently, PocketRN, AliveCor, and Hippocratic.ai. Dr Turakhia is an employee and corporate officer of iRhythm Technologies Inc.

Figure 2. Symptom-Rhythm Correlation: Arrhythmias and Ectopic Beats



[†] Any Arrhythmia[†] includes episodes of AF ≥ 30 sec, SVT ≥ 90 bpm & ≥ 4 beats, VT ≥ 100 bpm & ≥ 4 beats, any VF, Pause ≥ 3 sec, and/or AVB (any 2nd Deg or CHB).
^{††} Detections^{††} denotes number of patients with detected rhythm out of the total study population (N = 1,100,337), Daily Symptoms group (N = 229,005), or Non-Daily Symptoms group (N = 871,332).
^{†††} SRC^{†††} reflects % patients with ≥ 1 rhythm-correlated button press during the wear period from among patients with detection of each rhythm. Reported for total population, Daily and Non-Daily groups.