

# Arrhythmias in Patients with Daily vs. Non-Daily Symptoms Undergoing Long-Term Continuous Patch ECG Monitoring

Anthony J. Battisti PhD<sup>1</sup>, Relana Pinkerton PhD<sup>1</sup>, Vladimir Fokin PhD<sup>1</sup>, Brent Wright RN, DrPHc<sup>1</sup>, Mintu P. Turakhia MD, MAS<sup>1,2</sup>

<sup>1</sup>iRhythm Technologies, San Francisco, CA; <sup>2</sup>Stanford University, Stanford, CA

## Background

- Patch-based long-term continuous ambulatory ECG monitoring (LTCM) of 14 days has shown greater arrhythmia yield compared to <48 hour Holter.
- However, Holter remains in use for patients who report daily symptoms based on clinical or payor preferences.


## Objective

- We compared arrhythmia yield in patients based on symptom frequency.

## Methods

- Retrospective cohort study of Zio<sup>®</sup> monitor or XT LTCM (iRhythm Technologies, San Francisco, CA) devices prescribed in the US between June 2023 and July 2024 (Figure 1).
- Devices included a patient-activated button to document symptoms.
- Inclusions: Patients ≥18 years; Devices worn for >7 to 14 days.
- Exclusions: 100% AF.
- Data were stratified by observed symptom frequency, defined as the number of button presses/day:
  - Daily Symptoms (≥ 1/day), or
  - Non-Daily Symptoms (< 1/day)
- ECG data was analyzed via a deep-learned AI algorithm and confirmed by certified cardiographic technicians.
- Mean time to first episode and % detection within 48 hrs. were determined for each arrhythmia type.
- Odds Ratios were calculated to compare Daily vs. Non-Daily frequency groups within the first 48 hours of wear time. Logistic regression models were used to calculate odds ratios adjusting for age and sex difference.

## Results



## Cohort Description

- Cohort size of 1,100,337 patients:
  - 21% Daily Symptoms
  - 79% Non-Daily Symptoms
- Patients with Daily Symptoms were more likely:
  - To be younger (mean age 50.9 vs. 64.1 yrs.)
  - To be female (66.8% vs. 52.6%)
  - To have lower arrhythmia yield (69.1% vs. 80.9%)

## Daily vs. Non-Daily Symptoms Findings

- Patients reporting Daily Symptoms had a longer time to first detected arrhythmia (2.71 ± 3.26 vs. 2.01 ± 2.86 days) and first actionable arrhythmia (4.52 ± 3.96 vs. 4.20 ± 3.93 days) than patients with Non-Daily Symptoms (Table 1).
- However, Daily Symptom patients experienced a shorter mean time to first detected AF, sustained SVT, VF, and Pause episodes.

- For all arrhythmia types, mean time to first detected episode was >48 hours, regardless of symptom frequency.
- Percent detection through 48 hours varied by rhythm type and symptom frequency.
- 36.0% of Daily and 40.1% of Non-Daily patients with Actionable Arrhythmias were detected within 48 hours.
- Percentage of total yield detected within 48 hours was lower for Daily (59.4%) vs. Non-Daily (70.1%) patients.
- Percent detection of SVT, VT, and AV Block within 48 hours was lower for the Daily vs. Non-Daily group.
- Percent detection of AF, sustained SVT, VF and Pause within 48 hours was higher for the Daily group.
- Differences in 48-hour detection between Daily and Non-Daily Symptom patients remained consistent in the age and sex adjusted logistic regression models.

## Limitations

- Symptom frequency was determined from the number of button presses during the monitoring period.
- Symptoms may have been present but not reported by button press in some cases.

**Table 1. Arrhythmias Detected at 48 hours by Symptom Frequency**

Rhythm	Total Population		Detection by Symptom Frequency					
	Total Wear Detection n = 1,100,337	48 Hour Detection	Daily Symptoms (≥1 per day) n = 229,005		Non-Daily <sup>†††</sup> Symptoms (<1 per day) n = 871,332		Odds Ratio* (95% CI)	P-value
			Detection within 48 hours (%)	Days to First Episode (Mean ± SD)	Detection within 48 hours (%)	Days to First Episode (Mean ± SD)		
<b>Total Arrhythmia Yield<sup>†</sup></b>	78.5% (863,393)	68.2% (589,046/863,393)	59.4% (87,946/148,061)	2.71 ± 3.26	70.1% (501,100/715,332)	2.01 ± 2.86	0.63 (0.62-0.63)	<0.0001
<b>Actionable Arrhythmias<sup>††</sup></b>	35.8% (393,892)	39.5% (155,447/393,892)	36.0% (20,710/57,519)	4.52 ± 3.96	40.1% (134,737/336,373)	4.20 ± 3.93	0.84 (0.83-0.86)	<0.0001
<b>AF</b>	7.2% (79,329)	62.8% (49,808/79,329)	66.5% (6,281/9,450)	2.37 ± 3.37	62.3% (43,522/69,879)	2.68 ± 3.56	1.20 (1.15-1.26)	<0.0001
<b>SVT (All)</b>	71.3% (784,824)	66.5% (521,513/784,824)	58.4% (77,999/133,549)	2.77 ± 3.27	68.1% (443,514/651,275)	2.15 ± 2.94	0.66 (0.65-0.67)	<0.0001
<b>SVT Sustained<sup>††</sup></b>	6.5% (71,275)	20.6% (14,683/71,275)	21.3% (2,144/10,054)	6.42 ± 3.98	20.5% (12,539/61,221)	6.63 ± 4.05	1.05 (1.00-1.11)	0.053
<b>VT</b>	25.5% (280,287)	30.7% (86,174/280,287)	27.7% (11,427/41,251)	5.22 ± 3.91	31.3% (74,747/239,036)	4.95 ± 3.95	0.84 (0.82-0.86)	<0.0001
<b>VF</b>	0.01% (116)	18.1% (21/116)	20.0% (5/25)	5.25 ± 3.56	17.6% (16/91)	6.20 ± 3.63	1.17 (0.38-3.59)	0.781
<b>Pause</b>	3.2% (35,032)	36.9% (12,291/35,032)	38.0% (1,649/4,335)	4.36 ± 3.81	36.7% (11,272/30,697)	4.61 ± 3.94	1.06 (0.99-1.13)	0.092
<b>AV Block (Any 2<sup>nd</sup> Degree or CHB)</b>	1.5% (16,715)	47.0% (7,856/16,715)	43.8% (966/2,208)	4.08 ± 3.98	47.5% (6,890/14,507)	3.82 ± 3.94	0.86 (0.79-0.94)	0.001

<sup>†</sup>Arrhythmia Yield defined as any episode of AF ≥30 sec, SVT ≥90 bpm & ≥4 beats, VT ≥100 bpm & ≥4 beats, any VF, Pause ≥3 sec, and/or AVB (any 2<sup>nd</sup> Deg or CHB).

<sup>††</sup>Actionable Arrhythmia defined as AF ≥30 sec, SVT ≥90 bpm & ≥30s, VT ≥100 bpm & ≥4 beats, any VF, Pause ≥3 sec, and/or AVB (any 2<sup>nd</sup> Deg or CHB).

<sup>†††</sup>Includes 32.7% (284,784/871,332) of patients who reported no symptoms.

\*Odds Ratios and P-values compare rates of detection within 48 hours for Daily vs. Non-Daily Symptoms groups.

## Conclusions

- In patients reporting daily symptoms with actionable arrhythmias detected through >7 days of LTCM, no actionable arrhythmias were detected during the first 48 hours in 64% of cases.
- In most patients, Holter (<48 hours) is inadequate for arrhythmia detection, as daily symptom reporting may be uncorrelated to arrhythmia frequency.

## Disclosures

- AJ Battisti, R Pinkerton, V Fokin, and B Wright are employees of and have received equity from iRhythm Technologies, Inc.
- Dr. Turakhia has received equity from iRhythm, Connect America, Evidently, PocketRN, AliveCor, and Hippocratic.ai. Dr Turakhia is an employee and corporate officer of iRhythm Technologies Inc.