

**ZIO**<sup>®</sup> BY iRHYTHM

# Zio<sup>®</sup> ECG Utilization Software (ZEUS) Instructions for Use

**Rx**  
ONLY



Depictions of screens or reports included in this manual are examples only.

Features available are dependent on your account settings.

# Contents

---

DESCRIPTION.....	2
INDICATIONS FOR USE.....	3
CONTRAINDICATIONS .....	4
WARNINGS .....	4
PRECAUTIONS .....	5
NOTICE OF PRIVACY PRACTICES.....	5
ZIOSUITE APPLICATIONS.....	5
LOGIN .....	6
DASHBOARD.....	7
ENROLLMENT .....	10
PATIENT DETAIL.....	11
REPORT ACCESS .....	13
REPORT INTERPRETATION.....	14
USER MANAGEMENT .....	16
REPORT TYPES .....	17
RHYTHM DETECTION.....	23
ASYMPTOMATIC ARRHYTHMIA DETECTION.....	24
HEART RATE CALCULATIONS.....	25
MEASUREMENT INTERVALS.....	26
ALGORITHMS .....	27
SECURITY .....	30
SYMBOLS.....	35
TROUBLESHOOTING.....	37

## DESCRIPTION

---

The ZioSuite® web portal and ZioSuite® mobile application are for healthcare professionals to manage and streamline clinical workflows associated with the Zio service. The Zio system consists of a long-term cardiac recording using a prescribed Zio device in combination with Zio ECG Utilization Software (ZEUS), a software system utilizing proprietary deep-learned algorithms for data analysis.

The ZioSuite software module within the ZEUS System offers capabilities to enroll patients into the Zio service, access clinical reports, manage/perform clinical report interpretation, and administer user settings/access.

# INDICATIONS FOR USE

---

The device is intended to capture, analyze and report symptomatic and asymptomatic cardiac events and continuous electrocardiogram information for long-term monitoring. After wear, ECG data from compatible monitoring devices is processed and analyzed by the ZEUS System. A final report is generated on the beat-to-beat information from the entire ECG recording. For the Zio AT service, the ZEUS System supports the capture and analysis of automatically-detected arrhythmia events, as well as the analysis of uploaded patient-triggered events.

The ZEUS System is indicated for use on patients 18 years or older who may be asymptomatic or who may suffer from transient symptoms such as palpitations, shortness of breath, dizziness, light-headedness, pre-syncope, syncope, fatigue, or anxiety and patients who are asymptomatic. The reports are provided for review by the intended user to render a diagnosis based on clinical judgment and experience. It is not intended for use on critical care patients.

## **Intended Patient Population**

18 years or older

## **Intended Use Environment**

Professional use only

## CONTRAINDICATIONS

---

- Do not use the ZEUS System for critical care patients.
- Do not use ZEUS for patients with symptomatic episodes where variations in cardiac performance could result in immediate danger to the patient or when real-time or in-patient monitoring should be prescribed.
- Do not use ZEUS for patients with known history of life threatening arrhythmias.
- QT interval measurements from the ZEUS System are not intended to replace measurements from a 12-lead ECG. The ZEUS System does not independently diagnose drug-induced QT interval changes or QT interval prolongation.

## WARNINGS

---

- The Zio AT device has a maximum limit of transmitting symptomatic (patient-triggered) and asymptomatic (auto-triggered) cardiac events.
  - Asymptomatic cardiac events are transmitted until the maximum limit of 500 transmissions is reached. When the transmission limit is reached, the asymptomatic cardiac events are no longer sent for review during wear.
  - Symptomatic cardiac events are transmitted until the maximum limit of 100 transmissions is reached. When the transmission limit is reached, the symptomatic cardiac events are no longer sent for review during wear.

If a transmission limit is reached, the Zio AT monitor continues to record ECG data, including symptomatic and asymptomatic cardiac events. The data will be fully analyzed and included in the end-of-wear report.

When the patient is approaching a maximum transmission limit, Customer Care contacts the prescribing physician's office and the patient to send the patient an additional Zio AT device.

- If a Zio AT device is activated before completing patient registration, notifications of clinically actionable arrhythmias will be delayed. To avoid delays, complete the patient registration prior to activating the Zio AT device. A completed patient registration is the prescription order for continuous ambulatory electrocardiogram (ECG) monitoring.

## **PRECAUTIONS**

---

- ZEUS System QT interval measurements are likely to underestimate the global QT measurement from a 12-lead ECG. The user should consider this when interpreting the ZEUS System QT measurements.
- AF/AFL burden estimates presented during wear in the Daily Report are algorithm generated estimates. However, AF burden presented in the Final Report is human reviewed.

## **NOTICE OF PRIVACY PRACTICES**

---



iRhythm is committed to upholding patient privacy and protecting personal information, in particular Protected Health Information (PHI) collected and processed in conjunction with our Zio Service. We commit to complying with all applicable privacy laws and allowing patients to exercise their rights via their doctor. Our full Notice of Privacy Practices, found at <https://www.irhythmtech.com/content/privacy> describes our privacy practices, our legal duties, and patients' rights concerning PHI.

## **ZIOSUITE APPLICATIONS**

---

Prior to accessing ZioSuite, a user account must be setup by your site administrator. The ZioSuite mobile application can be accessed by iPhone devices at the Apple App store and for Android devices using Google Play.

# LOGIN

ZioSuite Website	
<ol style="list-style-type: none"><li>1. Open a browser and go to <a href="http://www.ziosuite.com">www.ziosuite.com</a>.</li><li>2. Enter email address and then select "Continue".</li><li>3. Enter password and then select "Log In".</li></ol>	
ZioSuite Mobile Application	
<ol style="list-style-type: none"><li>1. Select the ZioSuite app icon on your mobile phone.</li><li>2. Enter email address and password, and then select "Log In".</li></ol>	



# DASHBOARD

## ZioSuite Website

Based on your user role, a targeted dashboard will be displayed. Below are the short list and actions available.



### Allied Health Professionals

- Register Patient
- Posted Final Reports (Download, Print, or Archive Report)
- Unregistered Monitors
- Zio AT Trigger Limit [Zio AT only]\*
- Active Patients
- Reports Pending Interpretation (Print or Download)
- Transmission Reports (Download, Print, or Archive Report) [Zio AT Only]

### Physicians

- Reports Pending Interpretation (Assign, Print, Download)
- Posted Final Reports (Download, Print, or Archive Report)
- My Interpretation History (Print or Download Report)
- Zio AT Trigger Limit [Zio AT only]\*
- Active Patients
- Transmission Reports (Download, Print, or Archive Report) [Zio AT Only]

\*Displays names of patients wearing the Zio AT patch who are approaching or have reached the trigger limit.



### Zio AT Trigger Limit [Zio AT only]

- **Reached** indicates the patient has reached one of the thresholds for maximum transmissions.

- **Approaching** indicates the patient is nearing a maximum threshold.
- **Additional monitor** provides information about an additional monitor that may be sent to this patient.
  - **Yes** indicates an additional monitor is being sent to the patient.
  - **Pending Decision Contact Customer Care** indicates the patient is approaching a maximum threshold and Customer Care will contact the prescribing physician's office. If you have questions, contact Customer Care.
  - **Declined** indicates the prescribing physician's office declined the additional monitor.

When an additional Zio AT device is activated, the patient is removed from the trigger limit dashboard.

See pages 20 to 22 for more details on trigger limits.

⚠ Zio AT Trigger Limit		
PATIENT	LIMIT	ADDITIONAL MONITOR
Robert Garcia	<b>Reached</b> 06/20/2023	<b>Yes</b> Monitor being sent
Alexander Herbert	<b>Approaching</b>	<b>Declined</b> Account declined monitor
Marie Bostrom	Approaching	<b>Pending Decision</b> Contact Customer Care
6 more >		

## ZioSuite Mobile Application

Located at the bottom of the application, dashboard allows the user to switch to one of the following screens:

- **Reports:** List of accessible final reports.
- **Interpret:** List of reports available to the logged in user pending interpretation.
- **Transmissions (Zio AT Only):** List of Transmission Reports.
- **Patients:** List of patient records accessible to the logged in user.
- **Settings:** Change Password, PIN/Face ID Setup, Call/Email Customer Support




# ENROLLMENT

## ZioSuite Website

1. Select “+ Register Patient” on the Dashboard. Or, select “Patients>+Register Patient” from the menu.
2. Select the account location.
3. For “Enrollment Type,” select “In-Clinic” or “Home” and select “Next.”
4. Provide the enrollment details.

Entries in all fields are required unless “(Optional)” is indicated.

5. To submit a completed registration, select “Complete Registration”.
6. To save entries for later completion, create a draft registration and enter only the “Last Name”, “First Name”, “Sex”, and “Date of Birth”, then select “Save As Draft”.
7. To exit registration without saving changes, select “Cancel”.



Start Registration

Location

Traditional

Enrollment Type

In-Clinic  Home

New Registration

This registration is for the Clinic of Traditional. [Edit](#) [Add New Patient Location](#)

Section	Fields
Patient Information	Last Name, First Name, Sex (Male selected), Date of Birth, MRN (Optional)
Patient Contact	Primary Phone Number, Email, Address (Street, City, State, Zip)
Provider	Primary Care Provider, Referring Provider, Referring Provider (Optional)
Billing	Billing Method, Billing Address (Street, City, State, Zip)
ZioSuite Details	Member ID, Date of Birth, Sex, Date of Birth

Emergency Contact

Cancel

Save As Draft

Complete Registration

# PATIENT DETAIL

## ZioSuite Website

1. To access the detailed view of a patient record, select the “Patients>Active Patient or Patients>All Patients” menu option.
2. Scroll through the patient list and select the name of the patient record of interest. Note: If unable to locate the patient, type the name into the search window and select the name.
3. The following information is provided in the Patient Detail screen:
  - Patient Demographics
  - Prescriptions
  - Prescription Clinical Reports

Note: Clinical Report can be viewed by selecting the link.

The screenshot displays the ZioSuite website interface. The top navigation bar includes 'Dashboard', 'Registration', 'Patients', 'Reports', and 'Admin'. The main content area is titled 'Active Patients (10/24)' and shows a list of patients with columns for Name, Status, Location, and Date. Below the list, there is a 'Search ID' section with a search bar and a 'Go to Patient/Prescription' button. The 'Prescription Information' section is visible, showing details for a patient named 'Sally, D'. The table below lists prescriptions with columns for Prescription Name, Dosage Form, Strength, Unit, Start Date, and Prescription Status. The first row shows 'MORPHINE' with a strength of '10000000' and a status of 'Filled'.

Prescription Name	Dosage Form	Strength	Unit	Start Date	Prescription Status
MORPHINE	Tablet	10000000	mg	10/10/2020	Filled

## ZioSuite Mobile Application

1. To access the detailed view of a patient record, select the Patients tab.
2. Scroll through the patient list and select the name of the patient record of interest.


Note: If unable to locate the patient, type the name into the search window and select the name.

3. The following information is provided in the Patient Detail screen:

- Patient Demographics
- Prescriptions

4. Click on a Prescription to access associate reports and Rx details.


Note: Clinical Report can be viewed by selecting in the Report Name link.



# REPORT ACCESS


## ZioSuite Website

1. To access a clinical report, select the “Reports>REPORT TYPE” menu option.
2. Scroll through the report list and select “View Report” of the report of interest.
3. Upon selection, the Clinical Report will be visible with options to download the Report.



## ZioSuite Mobile Application




1. To access a clinical report, select the Report tab.
2. Scroll through the report list and select tab of interest. Note: If unable to locate the patient, type the patient name into the search window and select the associated tab.
3. Upon selection, the Clinical Report will be visible with options to view the Report and the associated interpretation.



# REPORT INTERPRETATION

## ZioSuite Website



1. To interpret a clinical report, select the “Reports>Pending Interpretation” menu option.
2. Scroll through the report list and select “Interpret Report” of the report of interest.
3. Upon selection, the Clinical Report will be visible with options to indicate “Agree with Findings”, “Insert Preliminary Findings”, as well as to enter physician interpretation.
4. Select “Accept & Sign” to incorporate an electronically signed interpretation into the report PDF.








## ZioSuite Mobile Application

1. To interpret a clinical report, select the Interpret tab.
2. Scroll through the report list and select report of interest.
3. Upon selection, the Clinical Report will be visible with options to indicate “Agree with Findings”, “Insert Preliminary Findings”, as well as to enter physician interpretation.
4. Select “Accept & Sign” to incorporate an electronically signed interpretation into the report PDF.













# USER MANAGEMENT

ZioSuite Website																					
<h3 style="text-align: center;">Add User</h3> <ol style="list-style-type: none"> <li>To add a user, select the “Admin&gt;Add New User” menu option.</li> <li>Fill-in the user registration details and select the “Register” option.</li> </ol>	 <p>The screenshot shows the ZioSuite Admin dashboard. The 'Admin' menu is open, highlighting 'Add New User'. Below, the 'User Information' form is visible, including fields for Name, Email, Password, and Role, along with 'User Permissions' and 'End Notification Settings' sections.</p>																				
<h3 style="text-align: center;">User List</h3> <ol style="list-style-type: none"> <li>To view the list of users, select the “Admin&gt;User Management” menu option.</li> </ol>	 <p>The screenshot shows the ZioSuite Admin dashboard with the 'User Management' menu option highlighted. Below, a table lists users with columns for Name, Email, Role, Status, and Action.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Name</th> <th>Email</th> <th>Role</th> <th>Status</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>Sample, Sample</td> <td>sample@ziosuite.com</td> <td>Physician</td> <td>DISABLED</td> <td>Yes</td> </tr> <tr> <td>John Doe</td> <td>john.doe@ziosuite.com</td> <td>Physician</td> <td>ENABLED</td> <td>Yes</td> </tr> <tr> <td>Jane Smith</td> <td>jane.smith@ziosuite.com</td> <td>Physician</td> <td>ENABLED</td> <td>Yes</td> </tr> </tbody> </table>	Name	Email	Role	Status	Action	Sample, Sample	sample@ziosuite.com	Physician	DISABLED	Yes	John Doe	john.doe@ziosuite.com	Physician	ENABLED	Yes	Jane Smith	jane.smith@ziosuite.com	Physician	ENABLED	Yes
Name	Email	Role	Status	Action																	
Sample, Sample	sample@ziosuite.com	Physician	DISABLED	Yes																	
John Doe	john.doe@ziosuite.com	Physician	ENABLED	Yes																	
Jane Smith	jane.smith@ziosuite.com	Physician	ENABLED	Yes																	
<h3 style="text-align: center;">Deactivate/Activate User Account</h3> <ol style="list-style-type: none"> <li>Select the “Admin&gt;User Management” menu option to view the list of users.</li> <li>Select the user name of interest.</li> <li>Select “Activate” to activate an inactive user.</li> <li>Select “Deactivate” to make a user inactive and unable to access ZioSuite web and mobile application.</li> </ol>	 <p>The screenshot shows the 'User Information' form for a user. Below the form, there are two large buttons: 'Activate' and 'Deactivate'.</p>																				

# REPORT TYPES

## Final Report




The Final Report is a comprehensive, end-of-wear PDF report showing the findings from the recording. Preliminary findings are provided for clinician review.

<p><b>Patient Demographics</b></p> <ul style="list-style-type: none"> <li>• Patient Name</li> <li>• Date of Birth</li> <li>• Patient ID</li> <li>• Gender</li> <li>• Primary Indication</li> <li>• Prescribing Clinician</li> </ul> <p><i>If applicable</i></p> <ul style="list-style-type: none"> <li>• Managing Location</li> <li>• Referring Physician</li> </ul>	<div style="text-align: center;">  <b>Zio XT Final Report for Report, Patient #16</b> </div> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 25%;"><b>Date of Birth</b> 12/12/67 (51 yrs)</td> <td style="width: 25%;"><b>Patient ID</b></td> <td style="width: 25%;"><b>Gender</b> Female</td> <td style="width: 25%;"><b>Primary Indication (R94.311)</b> Abnormal electrocardiogram</td> </tr> <tr> <td><b>Prescribing Clinician</b> Dr. E. Physician</td> <td colspan="3"><b>Managing Location</b> San Francisco</td> </tr> </table>	<b>Date of Birth</b> 12/12/67 (51 yrs)	<b>Patient ID</b>	<b>Gender</b> Female	<b>Primary Indication (R94.311)</b> Abnormal electrocardiogram	<b>Prescribing Clinician</b> Dr. E. Physician	<b>Managing Location</b> San Francisco																		
<b>Date of Birth</b> 12/12/67 (51 yrs)	<b>Patient ID</b>	<b>Gender</b> Female	<b>Primary Indication (R94.311)</b> Abnormal electrocardiogram																						
<b>Prescribing Clinician</b> Dr. E. Physician	<b>Managing Location</b> San Francisco																								
<p><b>Enrollment Details</b></p> <ul style="list-style-type: none"> <li>• Enrollment Period</li> <li>• Analysis Time</li> </ul>	<table style="width: 100%;"> <tr> <td style="width: 50%;"><b>Enrollment Period</b> 13 days 19 hours 03/22/19, 05:24am to 04/05/19, 12:40am</td> <td style="width: 50%;"><b>Analysis Time</b> 13 days 19 hours (after artifact removed)</td> </tr> </table>	<b>Enrollment Period</b> 13 days 19 hours 03/22/19, 05:24am to 04/05/19, 12:40am	<b>Analysis Time</b> 13 days 19 hours (after artifact removed)																						
<b>Enrollment Period</b> 13 days 19 hours 03/22/19, 05:24am to 04/05/19, 12:40am	<b>Analysis Time</b> 13 days 19 hours (after artifact removed)																								
<p><b>Arrhythmia Summary</b></p> <ul style="list-style-type: none"> <li>• Ventricular Fibrillation/Torsades de Pointes/Polymorphic VT</li> <li>• Ventricular Tachycardia</li> <li>• Pauses</li> <li>• Atrial Fibrillation/Flutter</li> <li>• AV Block (2<sup>nd</sup> Mobitz II, 3<sup>rd</sup>)</li> <li>• Supraventricular Tachycardia</li> </ul>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 80%;"><b>Ventricular Tachycardia (4 beats or more)</b></td> <td style="width: 20%;"><b>Episodes</b> 5</td> </tr> <tr> <td> <p>▼ Fastest VT (HR Range 135-150 bpm, Avg 142 bpm)</p>  </td> <td><b>HR Range</b> 116-150 bpm</td> </tr> <tr> <td></td> <td><b>Avg</b> 132 bpm</td> </tr> <tr> <td><b>Pauses (3 secs or longer)</b></td> <td><b>Episodes</b> 3</td> </tr> <tr> <td> <p>▼ Longest Pause (4.9 s, 12 bpm)</p>  </td> <td><b>Range</b> 3.9-4.9 s</td> </tr> <tr> <td><b>Atrial Fibrillation</b></td> <td><b>AF Burden</b> 37%</td> </tr> <tr> <td> <p>▼ Fastest AF (HR Range 61-154 bpm, Avg 100 bpm)</p>  </td> <td><b>Longest Duration</b> 1 d 19 h</td> </tr> <tr> <td></td> <td><b>HR Range</b> 50-154 bpm</td> </tr> <tr> <td></td> <td><b>Avg</b> 97 bpm</td> </tr> <tr> <td><b>AV Block (2<sup>nd</sup> Mobitz II, 3<sup>rd</sup>)</b></td> <td></td> </tr> <tr> <td style="height: 40px;">None found</td> <td></td> </tr> <tr> <td><b>Supraventricular Tachycardia (4 beats or more)</b></td> <td><b>None found</b></td> </tr> </table>	<b>Ventricular Tachycardia (4 beats or more)</b>	<b>Episodes</b> 5	<p>▼ Fastest VT (HR Range 135-150 bpm, Avg 142 bpm)</p> 	<b>HR Range</b> 116-150 bpm		<b>Avg</b> 132 bpm	<b>Pauses (3 secs or longer)</b>	<b>Episodes</b> 3	<p>▼ Longest Pause (4.9 s, 12 bpm)</p> 	<b>Range</b> 3.9-4.9 s	<b>Atrial Fibrillation</b>	<b>AF Burden</b> 37%	<p>▼ Fastest AF (HR Range 61-154 bpm, Avg 100 bpm)</p> 	<b>Longest Duration</b> 1 d 19 h		<b>HR Range</b> 50-154 bpm		<b>Avg</b> 97 bpm	<b>AV Block (2<sup>nd</sup> Mobitz II, 3<sup>rd</sup>)</b>		None found		<b>Supraventricular Tachycardia (4 beats or more)</b>	<b>None found</b>
<b>Ventricular Tachycardia (4 beats or more)</b>	<b>Episodes</b> 5																								
<p>▼ Fastest VT (HR Range 135-150 bpm, Avg 142 bpm)</p> 	<b>HR Range</b> 116-150 bpm																								
	<b>Avg</b> 132 bpm																								
<b>Pauses (3 secs or longer)</b>	<b>Episodes</b> 3																								
<p>▼ Longest Pause (4.9 s, 12 bpm)</p> 	<b>Range</b> 3.9-4.9 s																								
<b>Atrial Fibrillation</b>	<b>AF Burden</b> 37%																								
<p>▼ Fastest AF (HR Range 61-154 bpm, Avg 100 bpm)</p> 	<b>Longest Duration</b> 1 d 19 h																								
	<b>HR Range</b> 50-154 bpm																								
	<b>Avg</b> 97 bpm																								
<b>AV Block (2<sup>nd</sup> Mobitz II, 3<sup>rd</sup>)</b>																									
None found																									
<b>Supraventricular Tachycardia (4 beats or more)</b>	<b>None found</b>																								

<p><b>Heart Rate Summary</b></p> <ul style="list-style-type: none"> <li>Overall maximum, minimum and average heart rate</li> <li>Sinus maximum, minimum and average heart rate</li> </ul>	<p><b>Heart Rate</b></p> <table border="1"> <tr> <td><b>Overall</b></td> <td>Max</td> <td><b>154 bpm</b></td> <td>09:49am, 03/25</td> </tr> <tr> <td></td> <td>Min</td> <td><b>50 bpm</b></td> <td>11:59pm, 03/22</td> </tr> <tr> <td></td> <td>Avg</td> <td><b>78 bpm</b></td> <td></td> </tr> <tr> <td><b>Sinus</b></td> <td>Max</td> <td><b>96 bpm</b></td> <td>11:14am, 03/24</td> </tr> <tr> <td></td> <td>Min</td> <td><b>50 bpm</b></td> <td>11:59pm, 03/22</td> </tr> <tr> <td></td> <td>Avg</td> <td><b>66 bpm</b></td> <td></td> </tr> </table>	<b>Overall</b>	Max	<b>154 bpm</b>	09:49am, 03/25		Min	<b>50 bpm</b>	11:59pm, 03/22		Avg	<b>78 bpm</b>		<b>Sinus</b>	Max	<b>96 bpm</b>	11:14am, 03/24		Min	<b>50 bpm</b>	11:59pm, 03/22		Avg	<b>66 bpm</b>																					
<b>Overall</b>	Max	<b>154 bpm</b>	09:49am, 03/25																																										
	Min	<b>50 bpm</b>	11:59pm, 03/22																																										
	Avg	<b>78 bpm</b>																																											
<b>Sinus</b>	Max	<b>96 bpm</b>	11:14am, 03/24																																										
	Min	<b>50 bpm</b>	11:59pm, 03/22																																										
	Avg	<b>66 bpm</b>																																											
<p><b>Patient Events Summary</b></p> <ul style="list-style-type: none"> <li>Count of trigger events</li> <li>Count of diary events</li> <li>Range of findings for each event type</li> </ul>	<p><b>Patient Events</b></p> <p><b>Total Triggers: 2      Total Diaries: 1</b></p> <p>Findings within ± 45 sec of triggered events or diary entries:</p> <table border="1"> <thead> <tr> <th></th> <th>Range</th> <th>Trigger</th> <th>Diary</th> </tr> </thead> <tbody> <tr> <td>AF</td> <td>59-126 bpm</td> <td>✓</td> <td>✓</td> </tr> <tr> <td>Pause(s)</td> <td>3.9 s</td> <td>✓</td> <td></td> </tr> <tr> <td>Sinus</td> <td>56-73 bpm</td> <td>✓</td> <td></td> </tr> <tr> <td>SVE(s)</td> <td></td> <td>✓</td> <td></td> </tr> <tr> <td>VE(s)</td> <td></td> <td></td> <td>✓</td> </tr> </tbody> </table>		Range	Trigger	Diary	AF	59-126 bpm	✓	✓	Pause(s)	3.9 s	✓		Sinus	56-73 bpm	✓		SVE(s)		✓		VE(s)			✓																				
	Range	Trigger	Diary																																										
AF	59-126 bpm	✓	✓																																										
Pause(s)	3.9 s	✓																																											
Sinus	56-73 bpm	✓																																											
SVE(s)		✓																																											
VE(s)			✓																																										
<p><b>Ectopic Beats Summary</b></p> <ul style="list-style-type: none"> <li>Supraventricular Ectopy</li> <li>Ventricular Ectopy</li> </ul>	<p><b>Ectopics</b></p> <table border="1"> <thead> <tr> <th></th> <th>Rare &lt;1%</th> <th>Occasional 1% to ≤5%</th> <th>Frequent &gt;5%</th> </tr> </thead> <tbody> <tr> <td colspan="4"><b>Supraventricular Ectopy (SVE/PACs)</b></td> </tr> <tr> <td>Isolated</td> <td><b>Rare</b></td> <td>&lt;1.0%</td> <td>6723</td> </tr> <tr> <td>Couplet</td> <td><b>Rare</b></td> <td>&lt;1.0%</td> <td>141</td> </tr> <tr> <td>Triplet</td> <td><b>Rare</b></td> <td>&lt;1.0%</td> <td>9</td> </tr> <tr> <td colspan="4"><b>Ventricular Ectopy (VE/PVCs)</b></td> </tr> <tr> <td>Isolated</td> <td><b>Rare</b></td> <td>&lt;1.0%</td> <td>1716</td> </tr> <tr> <td>Couplet</td> <td><b>Rare</b></td> <td>&lt;1.0%</td> <td>192</td> </tr> <tr> <td>Triplet</td> <td><b>Rare</b></td> <td>&lt;1.0%</td> <td>26</td> </tr> <tr> <td>Longest Ventricular Bigeminy Episode</td> <td></td> <td></td> <td>0 s</td> </tr> <tr> <td>Longest Ventricular Trigeminy Episode</td> <td></td> <td></td> <td>0 s</td> </tr> </tbody> </table>		Rare <1%	Occasional 1% to ≤5%	Frequent >5%	<b>Supraventricular Ectopy (SVE/PACs)</b>				Isolated	<b>Rare</b>	<1.0%	6723	Couplet	<b>Rare</b>	<1.0%	141	Triplet	<b>Rare</b>	<1.0%	9	<b>Ventricular Ectopy (VE/PVCs)</b>				Isolated	<b>Rare</b>	<1.0%	1716	Couplet	<b>Rare</b>	<1.0%	192	Triplet	<b>Rare</b>	<1.0%	26	Longest Ventricular Bigeminy Episode			0 s	Longest Ventricular Trigeminy Episode			0 s
	Rare <1%	Occasional 1% to ≤5%	Frequent >5%																																										
<b>Supraventricular Ectopy (SVE/PACs)</b>																																													
Isolated	<b>Rare</b>	<1.0%	6723																																										
Couplet	<b>Rare</b>	<1.0%	141																																										
Triplet	<b>Rare</b>	<1.0%	9																																										
<b>Ventricular Ectopy (VE/PVCs)</b>																																													
Isolated	<b>Rare</b>	<1.0%	1716																																										
Couplet	<b>Rare</b>	<1.0%	192																																										
Triplet	<b>Rare</b>	<1.0%	26																																										
Longest Ventricular Bigeminy Episode			0 s																																										
Longest Ventricular Trigeminy Episode			0 s																																										
<p><b>Preliminary Findings</b></p> <ul style="list-style-type: none"> <li>Summary text of findings</li> </ul>	<p><b>Preliminary Findings</b></p> <p>Patient had a min HR of 50 bpm, max HR of 154 bpm, and avg HR of 78 bpm. Predominant underlying rhythm was Sinus Rhythm. 5 Ventricular Tachycardia runs occurred, the run with the fastest interval lasting 4 beats with a max rate of 150 bpm, the longest lasting 4 beats with an avg rate of 127 bpm. Episodes of Ventricular Tachycardia may be possible Atrial Fibrillation with aberrancy. Atrial Fibrillation occurred (37% burden), ranging from 50-154 bpm (avg of 97 bpm), the longest lasting 1 day 19 hours with an avg rate of 97 bpm. 3 Pauses occurred, the longest lasting 4.9 secs (12 bpm). Atrial Fibrillation and Pause were detected within +/- 45 seconds of symptomatic patient event(s). Isolated SVEs were rare (&lt;1.0%, 6723), SVE Couplets were rare (&lt;1.0%, 141), and SVE Triplets were rare (&lt;1.0%, 9). Isolated VEs were rare (&lt;1.0%, 1716), VE Couplets were rare (&lt;1.0%, 192), and VE Triplets were rare (&lt;1.0%, 26).</p>																																												
<p><b>Final Findings</b></p> <ul style="list-style-type: none"> <li>Clinician interpretation</li> </ul>	<p><b>Final Interpretation</b></p> <ol style="list-style-type: none"> <li>Agree with above interpretation</li> <li>Underlying Sinus rhythm with normal rates average =78/min</li> <li>5 runs of VT some of which could be AF with aberrancy</li> <li>Atrial fibrillation with 37% burden and longest run of 42 hours</li> <li>Pauses of up to 4.9 seconds likely post conversion related</li> <li>Triggered events consistent with AF, Pauses</li> </ol> <p>Electronically signed by Dr. Example Physician 04/12/19 06:18 PM (CT)</p> <p style="text-align: right;">SIGNATURE</p>																																												


# Transmission Report

Transmission Reports are provided for asymptomatic, symptomatic, scheduled, and baseline transmission events during the wear period of the Zio AT service only.

<p><b>Patient Demographics</b></p> <ul style="list-style-type: none"> <li>• Patient Name</li> <li>• Date of Birth</li> <li>• Patient ID</li> <li>• Gender</li> <li>• Primary Indication</li> <li>• Prescribing Clinician</li> </ul> <p><b>If applicable</b></p> <ul style="list-style-type: none"> <li>• Managing Location</li> <li>• Referring Physician</li> </ul>	<div style="text-align: right;">  <p><b>Zio AT Transmission Report for Report, Sample #4</b></p> </div> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 25%;">Date of Birth 12/12/67 (49 yrs)</td> <td style="width: 25%;">Patient ID</td> <td style="width: 15%;">Gender Female</td> <td style="width: 35%;">Primary Indication (I48.0) Paroxysmal atrial fibrillation</td> </tr> <tr> <td>Prescribing Clinician Dr. E. Physician</td> <td colspan="3">Managing Location INCC Lincolnshire</td> </tr> </table>	Date of Birth 12/12/67 (49 yrs)	Patient ID	Gender Female	Primary Indication (I48.0) Paroxysmal atrial fibrillation	Prescribing Clinician Dr. E. Physician	Managing Location INCC Lincolnshire		
Date of Birth 12/12/67 (49 yrs)	Patient ID	Gender Female	Primary Indication (I48.0) Paroxysmal atrial fibrillation						
Prescribing Clinician Dr. E. Physician	Managing Location INCC Lincolnshire								
<p><b>Enrollment Details</b></p> <ul style="list-style-type: none"> <li>• Enrollment Start</li> <li>• Date Received</li> <li>• Report Type</li> <li>• Transmission Count</li> </ul>	<table style="width: 100%;"> <tr> <td><b>Enrollment Start</b> 12/21/16, 01:34pm</td> <td><b>Date Received</b> 12/21/16, 09:46pm</td> </tr> <tr> <td><b>Report Type</b> Routine</td> <td><b>Transmission(s)</b> 1</td> </tr> </table>	<b>Enrollment Start</b> 12/21/16, 01:34pm	<b>Date Received</b> 12/21/16, 09:46pm	<b>Report Type</b> Routine	<b>Transmission(s)</b> 1				
<b>Enrollment Start</b> 12/21/16, 01:34pm	<b>Date Received</b> 12/21/16, 09:46pm								
<b>Report Type</b> Routine	<b>Transmission(s)</b> 1								
<p><b>Preliminary Findings</b></p> <ul style="list-style-type: none"> <li>• Summary text of findings</li> </ul>	<p><b>Preliminary Findings</b></p> <p>Atrial Fibrillation occurred continuously ranging from 152-187 bpm (avg of 170 bpm).</p>								
<p><b>Transmission Strip</b></p> <ul style="list-style-type: none"> <li>• 90-second transmission</li> </ul>	<div style="border: 1px solid #ccc; padding: 5px;"> <p><b>AF (152-187 bpm)</b> 12/21/16 09:29:09pm <span style="float: right;">= Patient Triggered Marker 02:56, 90s</span></p>  <p><small>Provided strip covers +5 seconds before and after the button areas.</small></p> </div>								
<p><b>Additional Strips</b></p> <ul style="list-style-type: none"> <li>• Enlarged segments of ECG</li> </ul>	<div style="border: 1px solid #ccc; padding: 5px;"> <p>12/21/16 09:28:58pm <b>Atrial Fibrillation (152-187 bpm)</b> <small>Associated with patient triggered event</small></p>  </div>								

# Daily Report

The Daily Report provides a recap of events detected and reported for the previous day. Daily Reports are available with the Zio AT service only.

<p><b>Patient Demographics</b></p> <ul style="list-style-type: none"><li>• Patient Name</li><li>• Date of Birth and Age</li><li>• Patient ID</li><li>• Gender</li><li>• Primary Indication</li><li>• Prescribing Clinician</li></ul> <p><b><i>If applicable</i></b></p> <ul style="list-style-type: none"><li>• Managing Location</li><li>• Referring Physician</li></ul>	<div data-bbox="473 386 904 506"><p><b>Zio AT Daily Report for Report, Sample #4</b></p><table border="1"><tr><td>Date of Birth 12/12/67 (49 yrs)</td><td>Patient ID</td><td>Gender Female</td><td>Primary Indication (I48.0) Paroxysmal atrial fibrillation</td></tr><tr><td>Prescribing Clinician Dr. E. Physician</td><td colspan="3">Managing Location INCC Lincolnshire</td></tr></table></div>	Date of Birth 12/12/67 (49 yrs)	Patient ID	Gender Female	Primary Indication (I48.0) Paroxysmal atrial fibrillation	Prescribing Clinician Dr. E. Physician	Managing Location INCC Lincolnshire						
Date of Birth 12/12/67 (49 yrs)	Patient ID	Gender Female	Primary Indication (I48.0) Paroxysmal atrial fibrillation										
Prescribing Clinician Dr. E. Physician	Managing Location INCC Lincolnshire												
<p><b>Enrollment Details</b></p> <p>The following information is included in the report header on the first page of the Daily Report:</p> <ul style="list-style-type: none"><li>• Enrollment Start</li><li>• Monitoring Date</li><li>• Number of days since activation</li><li>• Trigger limit indicators for the quantity of triggers utilized and available for each category:<ul style="list-style-type: none"><li>- Patient Trigger (symptomatic triggers)</li><li>- Auto Trigger (asymptomatic triggers)</li></ul></li></ul>	<div data-bbox="468 889 912 1019"><table border="1"><thead><tr><th>Enrollment Start</th><th>Monitoring Date</th><th>Day</th></tr></thead><tbody><tr><td>07/31/23, 08:00am</td><td>07/31/23</td><td><b>1</b></td></tr><tr><td><b>Patient Trigger</b></td><td>(16 of 100 used)*</td><td><input type="range"/></td></tr><tr><td><b>Auto Trigger</b></td><td>(0 of 500 used)*</td><td><input type="range"/></td></tr></tbody></table></div>	Enrollment Start	Monitoring Date	Day	07/31/23, 08:00am	07/31/23	<b>1</b>	<b>Patient Trigger</b>	(16 of 100 used)*	<input type="range"/>	<b>Auto Trigger</b>	(0 of 500 used)*	<input type="range"/>
Enrollment Start	Monitoring Date	Day											
07/31/23, 08:00am	07/31/23	<b>1</b>											
<b>Patient Trigger</b>	(16 of 100 used)*	<input type="range"/>											
<b>Auto Trigger</b>	(0 of 500 used)*	<input type="range"/>											

When the patient approaches a maximum trigger limit, the color of the bar changes and the report displays the following message for the category indicated by the alert icon (ⓘ):

Enrollment Start	Monitoring Date	Day
07/31/23, 08:00am	08/07/23	8
<b>Patient Trigger</b>	(70 of 100 used)*	
<b>Auto Trigger</b>	(0 of 500 used)*	

Your patient is approaching the max trigger limit for the category noted ⓘ. If the limit is reached, the monitor will continue to record but will not send event data for review during wear for this category.

When the patient reaches a maximum trigger limit, the report displays the following message for the category indicated by the caution icon (⚠):

<b>Patient Trigger</b>	(100 of 100 used)*	
<b>Auto Trigger</b>	(122 of 500 used)*	


Your patient has reached the max trigger limit for the category noted ⚠, which means the monitor will continue to record but will not send event data for review during wear for this category.

Note: The monitor continues to record data if the maximum limit is reached for a category. Event data will not be sent for review during wear for the category.

### Daily Summary

- Events across 24 hours
- Half-hourly AF burden estimates\*


Daily Summary for Thursday, January 27, 2022



## Wear Summary

- Events across the wear period
- Daily AF burden estimates\*
- Footnote explaining maximum trigger limits for the device

Wear Summary from Tuesday, January 18 to Thursday, January 27



Zio AT provides up to 14 days of continuous monitoring. The monitor has a maximum threshold to send data during wear for two categories of triggers: 100 Patient Triggers (button presses) and 500 Auto Triggers (auto detected ECGs). If the maximum trigger limit is reached for either category, the monitor will **NOT** continue to send event data for the category that reaches the maximum limit. The monitor will continue to send event data for the other category and the daily scheduled transmissions requested by the physician. The monitor will also continue to record ECG data, which will be analyzed and included in the end of wear Final Summary Report.

When a patient is approaching the limit for either category, iRhythm will send the patient an additional monitor and notify the account. Please instruct the patient to wear the replacement monitor. Refer to the Zio AT Clinical Reference Manual for more details.

## MD Notification Criteria

- A summary of the criteria which would prompt clinician notification

### MD Notification Criteria

Range	Duration	Notification	
AF	<40 bpm or >180 bpm	60s	During Business Hours Only
First Documentation of AF	Any	60s	During Business Hours Only
VT	>150 bpm	15s	Notify During Business Hours AND After Hours
SVT	>200 bpm	60s	Notify During Business Hours AND After Hours
Pauses	-	6s	Notify During Business Hours AND After Hours
AV Block (Chaf® Mobitz I)	Symptomatic	Any	Notify During Business Hours AND After Hours
AV Block (Chaf®)	Any	Any	Notify During Business Hours AND After Hours
Bradycardia	Symptomatic <40 bpm	30s	Notify During Business Hours AND After Hours

\*AF during periods of ventricular regularity or presence of Atrial Flutter may affect the estimated AF/AFL burden reported.



# RHYTHM DETECTION

---

The algorithm can analyze up to 14 days of data and detect the following rhythms:

- Pause  $\geq 3$  seconds
- Ventricular fibrillation
- Atrial fibrillation
- Complete heart block
- Second degree AV block type II
- Sinus rhythm (normal rhythm)
- Supraventricular tachycardia
- Ventricular bigeminy
- Ventricular tachycardia
- Ventricular trigeminy
- Second degree AV block type I
- Ectopic atrial rhythm
- Junctional rhythm
- Idioventricular rhythm

In addition to this list of arrhythmias, ECG findings in the end-of-wear report are observations from a qualified monitoring service technician.

# ASYMPTOMATIC ARRHYTHMIA DETECTION

Asymptomatic arrhythmia events, as detected and transmitted during the monitoring period, are defined by the following parameters:

Rhythm	Heart Rate	Duration
Atrial Fibrillation	≤40 bpm	≥60 seconds
	Between 40-180 bpm	≥60 seconds until first confirmation of AF
	≥180 bpm	≥60 seconds
Ventricular Tachycardia	≥120 bpm	≥30 seconds
	≥150 bpm	≥10 seconds
Supraventricular Tachycardia	≥180 bpm	≥60 seconds
Pause	-	≥4 seconds
	-	≥3 seconds back-to-back
Complete Heart Block	≤50 bpm	≥6 beats
Sinus Tachycardia	≥200 bpm	≥60 seconds
Sinus Bradycardia	≤30 bpm	≥60 seconds

For each of the arrhythmias listed above, the Zio AT patch will transmit up to four ECG (eight for CHB) strips per hour. Note: Following first confirmation of AF, an AF strip (any rate, ≥30 seconds) will be transmitted on a daily basis if detected on a subsequent day.

In addition to this list of arrhythmias, ECG findings in the reports provided during the monitoring period are observations from a qualified monitoring service technician.

Sensitivity (%) <sup>1</sup>		Positive Predictivity (%) <sup>1</sup>	
AHA	MIT-BIH	AHA	MIT-BIH
98.99	99.31	99.63	99.41

<sup>1</sup> TR01340.01 QRS detection (On file at iRhythm Technologies, Inc.)

# HEART RATE CALCULATIONS

---

<b>Episode Heart Rates</b>	<b>Max</b>	The maximum episode heart rate (i.e., maximum of all instantaneous heart rates within the episode)
	<b>Min</b>	The minimum episode heart rate (i.e., minimum of all instantaneous heart rates within the episode)
	<b>Avg</b>	The average episode heart rate (i.e., average of all instantaneous heart rates within the episode)
<b>Overall Rhythm Heart Rates</b>	<b>Max</b>	The maximum overall heart rate (i.e., maximum of all rhythm episode maximum heart rates within the record)
	<b>Min</b>	The minimum overall heart rate (i.e., minimum of all rhythm episode minimum heart rates exclusive of Pause heart rates within the record)
	<b>Avg</b>	The average overall heart rate (i.e., duration-weighted average of all rhythm episode heart rates within the record)

# MEASUREMENT INTERVALS

---

The Zio AT service may be configured to report QT interval measurements, along with a heart rate corrected QTc, from both symptomatic and asymptomatic transmissions using one of the following QT correction factors below.

QT Interval	The time period from the onset of the QRS interval waveform to the end/offset or the ensuing T-wave.
QTc Bazett	A correction factor for heart rate applied to QT as calculated by: $QT/\sqrt{RR}$
QTc Fridericia	A correction factor for heart rate applied to QT as calculated by: $QT/(RR)^{1/3}$
QTc Framingham	A correction factor for heart rate applied to QT as calculated by: $QT + 0.154(1-RR)$
QTc Hodges	A correction factor for heart rate applied to QT as calculated by: $QT + 1.75(\text{heart rate} - 60)$

QT Interval Measurements <sup>2</sup>			
Mean error	SD error	Min error	Max error
8.70 ms	7.25 ms	0.00 ms	45.00 ms

---

<sup>2</sup> TR01521.01 (On file at iRhythm Technologies, Inc.)

# ALGORITHMS

---

## ECG Deep Learning Analysis

The ECG Deep Learning Analysis algorithm (ECGDL) analyzes ECG recordings to provide beats, runs, rhythms, ECG segments, and heart rate detection.

## AutoTrigger Engine

The AutoTrigger Engine (ATE) is responsible for detecting asymptomatic arrhythmia events during the monitoring period.

## Algorithm Training

The source of training data for ECGDL and ATE algorithms are continuous cardiac recordings from compatible cardiac monitors. Training data is collected from thousands of recordings, which have already undergone Certified Cardiographic Technician (CCT) review.

## Algorithm Validation

Proprietary databases listed below were used to validate the algorithms. These consist of Zio ECG recordings from patients at least 18 years or older who may be asymptomatic or who may suffer from transient symptoms such as palpitations, shortness of breath, dizziness, light-headedness, pre-syncope, syncope, fatigue, or anxiety and patients who are asymptomatic.

## Bourn Database (ECGDL)

The Bourn Database, a proprietary, physician-validated database used for rhythm detection verification, has been developed from unique Zio ECG records for each rhythm class.

<b>Recording device</b>	Zio XT Patch, Zio SR Patch
<b>Channel(s)</b>	Single-lead ECG (modified lead II), chest electrodes with hydrogel
<b>Recording length</b>	5 minutes
<b>Environment</b>	Ambulatory
<b>Demographics</b>	<b>Age:</b> Median = 68 [25%,75%] = [55,77] <b>Gender:</b> 39% Female <b>Regional Demographics (USA):</b> <ul style="list-style-type: none"><li>• West: 44%</li><li>• Southwest: 5%</li><li>• Midwest: 18%</li><li>• Southeast: 22%</li><li>• Northeast: 11%</li></ul>

## ZBHT Database (ECGDL)

The ZBHT Database, a proprietary, physician-validated database, used for beat detection validation, has been developed from unique Zio ECG records for each beat class.

<b>Recording device</b>	Zio XT Patch, Zio AT Patch
<b>Channel(s)</b>	Single-lead ECG (modified lead II), chest electrodes with hydrogel
<b>Recording length</b>	30 seconds
<b>Environment</b>	Ambulatory
<b>Demographics</b>	<b>Age:</b> Median = 71 [25%,75%] = [60,78] <b>Gender:</b> 49% Female <b>Regional Demographics (USA):</b> <ul style="list-style-type: none"><li>• West: 36%</li><li>• Southwest: 6%</li><li>• Midwest: 19%</li><li>• Southeast: 24%</li><li>• Northeast: 15%</li></ul>

## QT Database (ECGDL)

The QT database, a proprietary, physician-validated database, used for QT measurement validation, has been developed from unique Zio ECG records.

<b>Recording device</b>	Zio XT Patch
<b>Channel(s)</b>	Single-lead ECG (modified lead II), chest electrodes with hydrogel
<b>Recording length</b>	90 seconds
<b>Environment</b>	Ambulatory
<b>Demographics</b>	<b>Age:</b> Median = 42 [25%,75%] = [31,62] <b>Gender:</b> 61% Female <b>Regional Demographics (USA):</b> <ul style="list-style-type: none"><li>• West: 34%</li><li>• Southwest: 2%</li><li>• Midwest: 23%</li><li>• Southeast: 25%</li><li>• Northeast: 16%</li></ul>

## ATFR Database (ATE)

The AT Full Recording (ATFR) database, used for asymptomatic event detection validation consists of multi-day ECG recordings obtained from the Zio AT Patch, along with CCT-reviewed reference labels obtained from the commercial Zio AT Service.

<b>Recording device</b>	Zio AT Patch
<b>Channel(s)</b>	Single-lead ECG (modified lead II), chest electrodes with hydrogel
<b>Recording length</b>	Up to 14 days (median: 13.8 days)
<b>Environment</b>	Ambulatory
<b>Demographics</b>	<b>Age:</b> Median=72 [25%,75%] = [60,80] <b>Gender:</b> 44.5% Female <b>Regional Demographics (USA):</b> <ul style="list-style-type: none"><li>• West: 18%</li><li>• Southwest: 6%</li><li>• Midwest: 16%</li><li>• Southeast: 27%</li><li>• Northeast: 33%</li></ul>

# SECURITY

---

iRhythm Technologies, Inc. uses industry best practices that ensure the confidentiality, integrity, and availability of data. Hosted at Amazon Web Services, our infrastructure is highly durable, scalable, and secure. We develop, manage, and maintain all proprietary software, systems, and associated security.

We are dedicated to exceeding our customer's expectations with respect to protected health information privacy and security by adhering to all relevant security requirements.

As participants in patient health care, we are committed to maintaining the privacy of Protected Health Information (PHI) as directed by applicable federal and state law. Our full Notice of Privacy Practices, found at [www.irhythmtech.com/content/privacy](http://www.irhythmtech.com/content/privacy) describes our privacy practices, our legal duties and rights concerning PHI.

## Certifications, Standards, Regulations



### SOC 2 Type II

Zio by iRhythm is SOC 2 Type II certified adhering to the AICPA's Trust Services Principles and Criteria for Security, Availability, Confidentiality and Privacy. The SOC 2 Type II is performed by an independent third-party and demonstrates iRhythm's commitment to Security and Privacy.



### HIPAA

The Health Insurance Portability and Accountability Act of 1996 (HIPAA) is a highly regulated and security-conscious statute in the healthcare industry. Zio by iRhythm is committed to maintaining HIPAA compliance and is regularly audited by independent third party assessors to help ensure we remain compliant.





**National Institute of Standards and Technology**

### **FIPS 140-2 Validation**

Zio by iRhythm has received the National Institute of Standards and Technology's (NIST) Federal Information Processing Standard (FIPS) 140-2 validation for data encryption. This achieves an added level of security required by specific government healthcare agencies and further demonstrates iRhythm's continued commitment to patient privacy and data security. Certificate number #3118.



### **CCPA**

The California Consumer Privacy Act (CCPA) is a state statute intended to enhance privacy rights and consumer protection for residents of California. Zio by iRhythm performs periodic independent third party Information Security / Data Privacy assessments to help with our compliance with requirements.



### **GDPR**

The General Data Protection Regulation (GDPR) is a regulation on data protection and privacy in the European Union. Zio by iRhythm receives regular independent third party assessments to help ensure we follow best practices in our efforts to comply with GDPR.



### **Privacy Shield**

Zio by iRhythm has chosen to continue our participation in the EU/US Privacy Shield Framework operated by the US Department of Commerce.

# Information Security

## Security

- Data encrypted in motion and at rest (HTTPS, AES-256)
- Role-based access controls
- 24/7 monitoring
- Regular penetration and vulnerability testing

## Cloud-Based

- AWS EC2 platform
- HL7-based EHR integration
- No on-premise hardware
- Highly scalable

## Availability

- Highly durable, geographically distributed architecture
- Scalable, virtualized server environment
- Redundant systems, no single point of failure
- Encrypted backups with offsite replication

## Auditing

- Comprehensive audit logging and alerting framework
- Activity tracking
- Regular risk assessments

## Policies and Procedures

- Extensive internal policy, procedure and operational controls
- Business Continuity Plan, including virtualization, cloud computing and dual site configuration
- Incident Response policy and procedures
- Business Associate Agreement with vendors that are involved with the delivery of the Zio Service.

## Cybersecurity Frequently Asked Questions

**1) What are the communication protocols used for ZioSuite website?**

ZioSuite runs on any current web browser, on any standard operating system. Information is encrypted and securely transmitted using industry-standard TLS 1.2 or greater protocol using regularly reviewed cypher suites.

**2) How should I maintain or update cybersecurity?**

ZioSuite is a Software as a Service (SaaS), cloud-based application with no on-premises hardware or software. ZioSuite.com runs on any current browser. It is recommend your browser is updated to the most current version before accessing ZioSuite.com.

**3) What actions should I take to ensure cybersecurity?**

It is recommended that users should only access ZioSuite.com using a trusted network connection. A trusted e-mail account should be used for the account creation process. Additionally, password best practices should be followed when creating an account. Passwords must be at least 8 characters; and should include at least 3 of (i) lower-case letters, (ii) upper-case letters, (iii) numbers, and (iv) special characters

**4) What types of cybersecurity events can be detected, and will I be notified?**





Data security is ensured using several security controls, including but not limited to, encryption at rest and in transit, firewalls, role-based access control, intrusion detection, audit logging and alerting, multi factor authentication, and complex password requirements. ZioSuite utilizes a comprehensive centralized logging and alerting system which records any data access, read, modification or removal, with associated user and timestamp. In addition, all account events are logged including, but not limited to, when and who created a customer account or related user accounts.


iRhythm complies with all relevant regulations for breach notification, particularly as applies to HIPAA and any security breach will be communicated with further instructions for next steps to mitigate.


**5) What should I do if a cybersecurity event is detected or suspected?**

If you were to suspect cybersecurity event, please contact Customer Support, available 24 hours a day 7 days a week, at 888-693-2401.

# SYMBOLS

Serial Number	 Catalogue number	<b>Rx ONLY</b> Prescription use only
 Manufacturer	 Date of manufacture	<b>QTY:</b> Net quantity of contents
 Consult instructions for use		

Symbol	Standard Reference	Standard Title	Symbol Title	Description/ Explanatory Text
	ISO 15223-1 Clause 5.1.1	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied	Manu- facturer	Indicates the medical device manufacturer
	ISO 7000-3082	Graphical symbols for use on equipment		
<b>Rx ONLY</b>	21 CFR 801.15(c) (1)(i)F	Labeling- Medical devices; prominence of required label statements	Prescription only	Requires prescription in the United States

<b>Symbol</b>	<b>Standard Reference</b>	<b>Standard Title</b>	<b>Symbol Title</b>	<b>Description/ Explanatory Text</b>
	ISO 15223-1 Clause 5.4.3	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 7000-1641	Graphical symbols for use on equipment		

# **TROUBLESHOOTING**

---

FOR CUSTOMER CARE, CALL **1.888.693.2401**




iRhythm Technologies, Inc., 699 8th St., Suite 600, San Francisco, CA 94103  
(888) 693-2401 | [www.irhythmtech.com](http://www.irhythmtech.com)

©2021-2025 iRhythm Technologies, Inc. All rights reserved.

Date of first issue: 2021-11 | Revised: 2025-02

LB10058.08

A decorative horizontal bar at the bottom of the page, transitioning from a light pink color on the left to a light blue color on the right.